



Scooter Lifts

Out-Sider® | Curb-Sider® | Joey® | Out-Rider™ | Chariot®

OUT-SIDER®

Lift capacity: 350 lb (159 kg)

SUVs, Vans, Pickups, Cars

3-Year Limited Warranty

Bruno's **Out-Sider** lets you easily transport your scooter or powerchair while retaining all your vehicle's seating and cargo space.

Lower the Out-Sider's self-leveling platform for easy drive-on access. Once your device is loaded, simply raise your mobility device into position where it will be locked in place using the Out-Sider's device-specific securement system.

The Out-Sider's platform automatically folds up to save space when unloaded, and it's low-profile design and adjustable license plate holder ensure maximum rear visibility for drivers.

Featuring an adaptable platform, the Out-Sider can be adjusted for a new mobility device or compatible vehicle in the future.

Options:

- Swing-Away for easy cargo access
- Large platform or center deck
- Low-profile third-brake light
- Securement belts
- Back-up camera
- Vinyl cover



Swing Away Option

JOEY™

Lift capacity: 350 lb (158 kg)
SUVs, Vans

3-Year Limited Warranty

The **Joey** offers a fully powered way to travel with your scooter or powerchair in any weather. Simply drive onto the self-leveling platform, press a button, and the Joey retracts into your vehicle's cargo area.

The Joey's compact design often preserves second-row seating, while a built-in exclusive safety barrier keeps the passenger area safe during sudden stops.

With illuminated controls for easy nighttime loading and the option to reinstall on another compatible vehicle, the Joey delivers convenience, safety, and long-lasting versatility.

■ Options:

- Wheel Chock Kit
(required for FWD/RWD powerchairs)
- Belt Kit



"I am able to take my scooter to the park, grocery shopping, church & even to restaurants."



James A.
Verified Bruno Customer





CURB-SIDER®

Lift capacity: 250 lb (114 kg), 400 lb (181 kg), 450 lb (203 kg)

SUVs, Vans, Pickups

3-Year Limited Warranty

The **Curb-Sider** leverages a versatile, compact design to effortlessly transport your mobility device while preserving seating and storage space.

Connect the Curb-Sider's docking device to your scooter or powerchair, press the "up" button, and the Curb-Sider lifts, rotates, and stores your mobility device for transport.

The Curb-Sider offers vehicle-specific bases for an optimized fit in your vehicle, as well as a telescoping head model that extends to clear large bumpers.

Adaptable and durable, the Curb-Sider can be reinstalled in another compatible vehicle or adjusted for a new mobility device.

■ Options:

- Quick release pins to fold down the lift head
- Exterior control (VSL-6000 only)
- Driver-side placement





OUT-RIDER®

Lift capacity: 100 lb (45 kg), 200 lb (90 kg), 350 lb (159 kg)

Pickup Box Size: 5½ ft, 6½ ft, 8 ft

3-Year Limited Warranty

The **Out-Rider** is designed to make lifting and transporting a mobility device in your truck easy.

Connect the Out-Rider's docking device to your scooter, powerchair, or wheelchair, press a button on the hand-held control and watch as the Out-Rider raises, rotates, and stows your mobility device with ease.

The Out-Rider meets or exceeds all applicable vehicle safety standards, and can be easily transferred to another compatible pickup for continued use.

■ Option:

- Passenger side model





CHARIOT®

Lift capacity: 350 LB (159 KG)

Sedans, SUVs, Vans

3-Year Limited Warranty

The **Chariot** connects to your hitch and delivers convenient transport of your scooter or powerchair outside your vehicle.

The Chariot stays perfectly in line with your vehicle while in drive or reverse thanks to 360-degree swiveling wheels.

Drive your scooter onto the self-leveling, traction-coated platform and press a button to raise or lower the lift. Self-tensioning retractable securement belts help keep your mobility device stable during travel.


The illuminated license plate holder, brake lights, running lights, and turn signals increase visibility. The Chariot's adjustable platform fits most mobility devices and can be reinstalled on other compatible vehicles.

Options:

- Platform corner markers
- Back-up camera
- Vinyl cover
- Spare tire



Scan for more
info & videos



Bruno's mission is to help make life easier for people with mobility challenges by delivering the highest quality stair lifts, scooter lifts, platform lifts and home elevators.