



ELAN

Straight Indoor Stair Lift

 **BRUNO****MADE IN USA**



BRUNO ELAN

#1 Selling Stairlift Made in the USA

Easy to Use

The Elan smoothly glides up and down your staircase with the press of a single button.


Space Saving

The Elan's purpose-built vertical rail maximizes open space on steps for family, friends and pets. The armrests, seat and footrest fold up, creating a compact storage profile.

Dependable

The Elan is battery-powered and continuously charges anywhere it stops along the rail, delivering reliable operation even in the event of a power outage.





AVAILABLE OPTIONS

Make rotating the Elan seat at the top of the steps automatic with the selection of the **power swivel** option.

Easily fold the footrest up into the storage position by choosing the popular **power footrest** option.

Maintain open space at the bottom of the staircase with the addition of a manual or **power folding rail** option.

Not all options shown.





TECH SPECS

Lift Capacity:

300 lb (136 kg)

Seat Size:

14.1 in (358 mm) deep x
15.38 in (391 mm) wide

Width Between Armrests:

20.5 in (520 mm)

Footrest Size:

10.5 in long x 10.75 in wide
(267 mm x 276 mm)

Vertical Rail:

Installs within 5 in (127 mm) of the wall

WARRANTY

Limited Lifetime Warranty on Major

Components (motor, gearbox, drivetrain and stationary rail system) for the lifetime of the Bruno stair lift as long as the original purchaser owns the stair lift.

2-Year Warranty on all parts after date of installation by an Authorized Bruno Dealer.

30-Days Labor Coverage after date of installation by an Authorized Bruno Dealer for any manufacturing defect.



Scan for more
info & videos



Bruno's mission is to help make life easier for people with mobility challenges by delivering the highest quality stair lifts, scooter lifts, platform lifts and home elevators.

© Bruno Independent Living Aids, Inc. 2026
2026-05-05-elan-mkt-sl-6.indd