



ELITE

Curved Indoor Stair Lift





ACCESS ALL LEVELS OF **YOUR HOME**

Experience the ultra-smooth ride of a Bruno Elite curved stair lift. Made in USA and handcrafted for a precise fit, the Elite's custom rail maximizes space and blends seamlessly with your staircase.

Bruno's flagship stair lift, the Elite is built for lasting independence at home - featuring high-performance materials, heavy-duty capacity and a Limited Lifetime Warranty.



NORTH AMERICA'S PREMIER CURVED STAIR LIFT

Handcrafted for Your Home

Elegant, custom-made rail delivers a precision fit and smooth ride while the covered gear rack provides a clean, graceful appearance.

Easy to Use

Armrest control glides the stair lift up and down. Includes two wireless remotes.

Customizable

Change seat height, footrest height and armrest width for optimal fit.

Saves Space

Vertical rail maximizes open space on steps. Armrests, seat and footrest fold flat when not in use.

Dependable

Battery-powered design delivers top-to-bottom operation even in the event of a power outage.

Peace of Mind

Obstruction sensors and seatbelt enhance safety during travel. Swiveling seat locks in place for safe exit away from stairs.



MADE IN USA

Bruno stair lifts are engineered and manufactured in Wisconsin. Each Elite curved rail is handcrafted to ensure an optimal fit to your staircase and exceptional ride quality.



"I love my Bruno stair lift. It works so smoothly. **It is truly a blessing.**"



Sharon W.
Bruno Stair Lift Customer

AVAILABLE OPTIONS

Select from four additional upholstery color options to blend your stair lift perfectly into your home's decor.

Add extended rail positions to park the stair lift away from the steps or out of view, or an optional mid-park charging station for stairs with middle landings.

Eliminate the need to lean over with a power folding footrest. Simply fold the seat and the footrest will automatically fold, too.

Choose the power swivel option to automatically rotate your seat with the push of a button.

Opt for a wider seat or larger footrest to further customize your Elite stair lift chair.

Not all options shown.





More Info
& Video

TECH SPECS

Lift Capacity:

400 lb (181 kg)

Standard Seat Size:

18 in wide x 15 in deep (457 mm x 381 mm)

Larger Seat Option:

22.5 in wide x 16.5 in deep
(572 mm) x (429 mm)

Standard Footrest Size:

10.5 in long x 9.4 in wide
(260 mm x 232 mm)

Larger Footrest Option:

12 in long x 11 in wide (305 mm x 279 mm)

Curved Rail:

Installs as close as 7.5 in (191 mm) to the wall or edge of staircase

WARRANTY

Limited Lifetime Warranty on Major

Components (motor, gearbox, drivetrain and stationary rail system) for the lifetime of the Bruno stair lift as long as the original purchaser owns the stair lift.

2-Year Warranty on all parts after date of installation by an Authorized Bruno Dealer.

30-Days Labor Coverage after date of installation by an Authorized Bruno Dealer for any manufacturing defect.

ABOUT BRUNO

Veteran-Founded | Family-Owned | Made in USA

Bruno's mission is to help make life easier for people with mobility challenges by delivering the highest quality stair lifts, scooter lifts, platform lifts and home elevators.



© Bruno Independent Living Aids, Inc. 2025
2025-09-09-elite-curve-mkt-sl3.indd