



ELITE

Straight Indoor Stair Lift



MADE IN USA



BRING HOME **COMFORT & ACCESS**

Enjoy modern styling and premium performance with Bruno's top-of-the-line Elite straight stair lift.

The Elite delivers a smooth ride from start to finish featuring high-quality materials and Made in USA quality. With several customizable options, the Elite is ready to meet your needs and blend seamlessly into your home.

BRUNO'S PREMIUM STRAIGHT STAIR LIFT

Personalized Fit

Generous padding, available power options, seat color upgrades and adjustable armrest, seat and footrest settings combine to deliver unmatched comfort and style.

Easy to Use

Armrest control glides the stair lift up and down. Includes two wireless remotes.

Saves Space

Vertical rail with hidden gear rack maximizes open space on steps and delivers a clean appearance. Armrests, seat and footrest fold flat when not in use.

Dependable

Continuously charging, battery-powered design delivers top-to-bottom operation even in the event of a power outage.

Peace of Mind

Obstruction sensors and seatbelt enhance safety during travel. Swiveling seat locks in place for safe exit away from stairs.

Limited Lifetime Warranty

The Bruno Elite is designed for lasting durability. Enjoy a Limited Lifetime Warranty on your new Elite for as long as you own it.



"I love my Bruno stair lift. It works so smoothly. **It is truly a blessing.**"



Sharon W.
Bruno Stair Lift Customer





AVAILABLE OPTIONS

Select from four additional upholstery color options to blend your stair lift perfectly into your home's decor.

Create more open space at the bottom of your staircase by adding a manual or automatic folding rail option.

Eliminate the need to lean over with a power folding footrest. Simply fold the seat and the footrest will automatically fold, too.

Choose the power swivel option to automatically rotate your seat with the push of a button.

Opt for a wider seat or larger footrest to further customize your Elite stair lift chair.

Not all options shown.





More Info
& Video

TECH SPECS

Lift Capacity:

400 lb (181 kg)

Standard Seat:

18 in wide x 15 in deep (457 mm x 381 mm)

Larger Seat Option:

22.5 in wide x 16.5 in deep
(572 mm) x (429 mm)

Standard Footrest:

10.5 in long x 9.4 in wide
(260 mm x 232 mm)

Larger Footrest Option:

12 in long x 11 in wide (305 mm x 279 mm)

Vertical Rail:

Installs within 6 in (152 mm) of the wall

WARRANTY

Limited Lifetime Warranty on Major

Components (motor, gearbox, drivetrain and stationary rail system) for the lifetime of the Bruno stair lift as long as the original purchaser owns the stair lift.

2-Year Warranty on all parts after date of installation by an Authorized Bruno Dealer.

30-Days Labor Coverage after date of installation by an Authorized Bruno Dealer for any manufacturing defect.

ABOUT BRUNO

Veteran-Founded | Family-Owned | Made in USA

Bruno's mission is to help make life easier for people with mobility challenges by delivering the highest quality stair lifts, scooter lifts, platform lifts and home elevators.



© Bruno Independent Living Aids, Inc. 2025
2025-09-09-elite-indoor-straight-mkt-sl-2.indd